

Pride® Lifts

Mobility Made Easy.

Outlander Full Platform

The **Outlander Full Platform** Exterior Lift for scooters and power chairs is the most universal platform lift available. The simple hitch-based installation requires little to no modification to your vehicle and allows for easy removal. The easy drive-on/drive-off platform enables loading from either end of the lift and takes the hassle out of loading your mobility device in a busy parking area. Consumers have the choice of a one-hand retractable strap kit, scooter lock-down arm or power chair hold-down bar.

BASE WEIGHT ¹	82 lbs.
MAXIMUM SCOOTER/PC OVERALL WEIGHT	350 lbs. ²
MAXIMUM PRODUCT WIDTH	28"
OVERALL PLATFORM LENGTH	43", 42" usable
OVERALL PLATFORM WIDTH	30"
MAXIMUM TRAVEL DISTANCE FROM THE GROUND	21"
MOTOR	Fully sealed 12 Volt DC
WARRANTY	3-Year Transferable
HITCH CLASS ³	II, III, IV
OPTIONS	<ul style="list-style-type: none">• Battery pack w/ charger• Swing-away adaptor• Scooter hold-down arm• Power chair hold-down bar• Strap kit
VEHICLE PLATFORM	Full-size van, Minivan, Truck, Sedan, SUV and Crossover
LIFTS CURRENT PRODUCTS ^{1,2}	All Pride® Scooters & Pride®/Quantum® PWCs with contoured, Specialty & Synergy® seating And equivalent competitor products Except: Pursuit®, Pursuit® XL, Pursuit® Sport 36 Volt, Sport Rider, Maxima 3 & 4, 1420, 1450 & 1650

1. Vehicle dependent.
2. Lifting weight does not include accessories, upgraded seating and/or hardware.
3. Hitch weight carrying capacity varies by vehicle. Call for details.
4. Does not include any accessories or hardware.

To find a compatible lift for your vehicle and mobility product, use our new

LIFT FINDER

pridemobility.com/lift-finder

FEATURES AND BENEFITS

- 350 lb. lifting capacity
- Compatible with Class II, III or IV hitches
- LED puddle light for ease of use at night
- Toggle switch for ease of up and down use
- Integrated license plate bracket
- Phosphate undercoating for weather resistance
- LED indicator for On/Off
- Manual override with tool kit
- Adjustable wheel chocks (tire stops that can be adjusted to wheelbase of PMD)
- Standard hardware kit



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.