



Jazzy® Elite ES

Jazzy® Elite ES-1

Jazzy® Elite ES Portable

	Jazzy® Elite ES	Jazzy® Elite ES-1	Jazzy® Elite ES Portable
Code coverage	K0822/K0823	K0816	K0820/K0821
Weight capacity	300 lbs.	300 lbs.	300 lbs.
Maximum speed	Up to 4 mph	Up to 4 mph	Up to 4 mph
Ground clearance	1.5"	0.75"	1.5"
Turning radius	24.75"	23.5"	25.5"
Length	34.75"	34.75"	34.75"
Width	22.75"	22.75"	22.75"
Drivetrain	Two-motor, in-line, front-wheel drive	Two-motor, in-line, front-wheel drive	Two-motor, in-line, front-wheel drive
Standard electronics	40A, PG GC3	40A, PG GC3	40A, PG GC3
Battery charger	Off-board, 3 amp	Off-board, 3 amp	Off-board, 3 amp
Per-charge range	Up to 15 miles	Up to 12 miles	Up to 11 miles
(2) Batteries	12 volt, U-1	12 volt, U-1	12 volt, U-1
Battery weight	24.5 lbs. each	18 lbs. each	18.4 lbs. each
Base weight	68.5 lbs.	68.5 lbs.	68.6 lbs.
Standard seating	18" x 18" comfort high-back seat with headrest	18" x 18" comfort high-back seat with headrest	18" x 18" comfort high-back seat with headrest
Disassembly	No	No	Yes

Options

- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

Colors



Jazzy Red



Jazzy Blue



Your Group 1 and 2 Power Wheelchair Solutions



Jazzy[®]

ELITE ES SERIES

Jazzy[®] Elite ES

Jazzy[®] Elite ES-1

Jazzy[®] Elite ES Portable



Your Group 1 and 2 Power Wheelchair Solutions

The Jazzy[®] Elite ES Series delivers a reliable blend of power, performance and style. In-line, front-wheel drive technology gives the Elite ES Series excellent stability and maneuverability for solid performance indoors and out.



Jazzy[®]

ELITE ES

Features

- In-line, front-wheel drive technology
- 40 amp, PG GC 3 controller
- Simple main frame design for easy serviceability
- 18" x 18" comfort high-back seat with headrest (standard)
- Jazzy[®] Armor[®] controller guard to protect from daily wear and tear



Jazzy[®]

ELITE ES-1

Features

- In-line, front-wheel drive technology
- 40 amp, PG GC 3 controller
- Simple main frame design for easy serviceability
- 18" x 18" comfort high-back seat with headrest (standard)
- Jazzy[®] Armor[®] controller guard to protect from daily wear and tear



Jazzy[®]

ELITE ES PORTABLE

Features

- Easy disassembly for convenient transport
- In-line, front-wheel drive technology
- 40 amp, PG GC 3 controller
- Simple main frame design for easy serviceability
- 18" x 18" comfort high-back seat with headrest (standard)

