

Quantum® Q6 Edge 2.0 with
Static Contoured Seating
Group 3



QUANTUM REHAB®



Discount: _____
Terms: _____
Freight: _____
Taken By: _____

Account Number: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
Email: _____
Date: _____
PO#: _____
Marked For: _____
Ship to Address: _____
City: _____ State: _____ Zip: _____

Note Section

For Quantum Use:

For Provider Use:


Set-to-Spec

This section must be filled out, otherwise the selection will default to the mid-range measurements.


Provide legrest height.

Set Legrest Length Left: _____ Right: _____

Height: _____ **Weight:** _____



- A. Shoulder Width _____
- B. Chest Width _____
- C. Hip Width _____
- D. External Knee Width _____



- E. Maximum Sitting Height _____
- F. Shoulder Height _____
- G. Axilla Height _____
- H. PSIS Height _____
- I. Trunk Depth _____
- J. Ischial Depth _____
- K. Thigh Depth _____
- L. Lower Leg Length _____
- M. Elbow Height _____

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.



Quantum Q6 Edge 2.0 with Contoured Seating

Description	Price	HCPCS Code
Base Model		
<input type="checkbox"/> Quantum Q6 Edge 2.0 3S-SS <i>300 lbs. wt capacity. Comes standard with a solid seat pan, 6 mph programmable motors, 8 amp off board charger. Dimensions of base without legrests are 24"W x 35.5"L</i>	\$7,595	K0848
<input type="checkbox"/> Quantum Q6 Edge 2.0 3S-C <i>300 lbs. wt capacity. Comes standard with a captain's seat (contoured seating), 6 mph programmable motors, 8 amp off board charger. Dimensions of base without legrests are 24"W x 35.5"L</i>	\$7,595	K0849
<input type="checkbox"/> Quantum Q6 Edge 2.0 3SP-SS <i>For use with PES. 300 lbs. wt capacity. Comes standard with a solid seat pan, 6 mph programmable motors, 8 amp off board charger. Dimensions of base without legrests are 24"W x 35.5"L</i>	\$7,595	K0856

Colors					
Metallic Colors:					
<input type="checkbox"/> Red Corvette	<input type="checkbox"/> Back in Black	<input type="checkbox"/> Pretty in Pink	<input type="checkbox"/> Lime Light	<input type="checkbox"/> Silver Bullet	<input type="checkbox"/> Orange Crush
<input type="checkbox"/> Electric Blue	<input type="checkbox"/> White Satin	<input type="checkbox"/> Deep Purple	<input type="checkbox"/> Grey Street	<input type="checkbox"/> Midnight Blue	<input type="checkbox"/> Yellow Blaze

Casters	
<input type="checkbox"/> 6" Solid Casters	Standard
<input type="checkbox"/> 8" Solid Casters	\$150




Base Mounted Foot Platform	
<input type="checkbox"/> Base Mounted Foot Platform	Standard
<input type="checkbox"/> Bell Shaped Foot Platform <i>(10.12" Tapered to 14"W x 13"D)</i>	\$235

Drive Wheels					
<input type="checkbox"/> 3" x 14" Flat Free Black Tires	Standard				
<input type="checkbox"/> 3" x 14" Pneumatic Black Tires	No Charge				
Rim Color:					
<input type="checkbox"/> Silver <i>(Standard)</i>	<input type="checkbox"/> Matching <input type="checkbox"/> Non-Matching				
Rim Color: \$50 <i>(Available through SOF, not available with 4"W tires)</i>					
<input type="checkbox"/> Red Corvette	<input type="checkbox"/> Back in Black	<input type="checkbox"/> Pretty in Pink	<input type="checkbox"/> Lime Light	<input type="checkbox"/> Midnight Blue	<input type="checkbox"/> Yellow Blaze
<input type="checkbox"/> Electric Blue	<input type="checkbox"/> White Satin	<input type="checkbox"/> Deep Purple	<input type="checkbox"/> Grey Street	<input type="checkbox"/> Orange Crush	

Batteries		
<input type="checkbox"/> NF-22, 50 Amp Batteries		
<input type="checkbox"/> Maintenance Free AGM <i>(UPG)</i>	\$295 ea.	E2361
<input type="checkbox"/> Maintenance Free AGM <i>(MK)</i>	\$295 ea.	E2361
<input type="checkbox"/> Maintenance Free Gel Cell <i>(UPG)</i>	\$295 ea.	E2361
<input type="checkbox"/> Maintenance Free Gel Cell <i>(MK)</i>	\$295 ea.	E2361
<input type="checkbox"/> No Batteries		
Installation:		
<input type="checkbox"/> Batteries Installed	<input type="checkbox"/> Do Not Install	







LED Light Package	
<input type="checkbox"/> Through Q-Logic 2 EX or NE+ 8-Key w/ Lights	\$995

Fender Lights	
<input type="checkbox"/> Fender Lights through Single Switch Only	\$275
<input type="checkbox"/> Fender Lights through Controller/Single Switch Combination	\$250

Non-Expandable Controls	
<input type="checkbox"/>  75 Amp NE	Standard
<input type="checkbox"/>  75 Amp NE+ 6-Key <i>For up to 2 Actuators through joystick</i>	\$200
<input type="checkbox"/>  75 Amp NE+ 8-Key <i>For Lights</i>	\$400



Quantum Q6 Edge 2.0 with Contoured Seating

Description	Price	HCPCS Code
Expandable Controls		
<input type="checkbox"/>  100 Amp Q-Logic 2 EX Controller <i>Required with Q-Logic Electronics</i>	\$695	
<input type="checkbox"/> Expandable Controller Harness <i>Required with Q-Logic Electronics</i>	\$425	
<input type="checkbox"/>  Q-Logic 2 EX Joystick	\$700	
<input type="checkbox"/>  Q-Logic 2 EX Joystick <i>For Lights</i>	\$700	
Q-Logic Color Shroud Options:		
<input type="checkbox"/> ● Black Soft Touch <i>Standard</i>	<input type="checkbox"/> Matching	<input type="checkbox"/> Non-Matching
Joystick Color Shroud: \$50		
<input type="checkbox"/> ● Red Corvette	<input type="checkbox"/> ● Back in Black	<input type="checkbox"/> ● Pretty in Pink
<input type="checkbox"/> ● Lime Light	<input type="checkbox"/> ● Silver Bullet	<input type="checkbox"/> ● Orange Crush
<input type="checkbox"/> ● Electric Blue	<input type="checkbox"/> ● White Satin	<input type="checkbox"/> ● Deep Purple
<input type="checkbox"/> ● Grey Street	<input type="checkbox"/> ● Midnight Blue	<input type="checkbox"/> ● Yellow Blaze
<input type="checkbox"/>  Stand Alone Joystick with External Charger Port	\$700	
<i>Features:</i>		
<ul style="list-style-type: none"> • Up to 3 drive profiles and 1 seating profile when used alone. • Use in conjunction w/ enhanced display for up to 5 profiles. • Auxiliary functions including Bluetooth Mouse and IR. 		
<input type="checkbox"/>  Q-Logic 2 Enhanced Display Kit	\$1,884.93	
<i>Includes:</i>		
<ul style="list-style-type: none"> • Enhanced Display • Flexible Gooseneck • Gooseneck Mounting Bracket 		
<i>Features:</i>		
<ul style="list-style-type: none"> • Specialty controls interface w/ 3.5" LCD display • Up to 5 profiles to control drive, seat and auxiliary functions • Integrated single switch scanner, Bluetooth mouse emulator and IR controls • Bluetooth connectivity for mobile devices 		
Enhanced Display Color Shroud Options:		
<input type="checkbox"/> ● Black SoftTouch <i>Standard</i>	<input type="checkbox"/> Matching	<input type="checkbox"/> Non-Matching
Display Color Shroud: \$50		
<input type="checkbox"/> ● Red Corvette	<input type="checkbox"/> ● Back in Black	<input type="checkbox"/> ● Pretty in Pink
<input type="checkbox"/> ● Lime Light	<input type="checkbox"/> ● Silver Bullet	<input type="checkbox"/> ● Orange Crush
<input type="checkbox"/> ● Electric Blue	<input type="checkbox"/> ● White Satin	<input type="checkbox"/> ● Deep Purple
<input type="checkbox"/> ● Grey Street	<input type="checkbox"/> ● Midnight Blue	<input type="checkbox"/> ● Yellow Blaze
<input type="checkbox"/>  USB Mobile Device Charger	\$175	
<i>Available with NE, NE+ and Q-Logic 2 EX controllers; Can only charge 1 device at a time. A multiplier harness is required with all Static and Tilt only seating systems. The selection of more than 3 modules requires an additional multiplier.</i>		

Joystick Handles		
<input type="checkbox"/> Standard Carrot Handle	Standard	
<input type="checkbox"/> Mushroom Handle	\$71	E2323
<input type="checkbox"/> T-Bar	\$107	E2323
<input type="checkbox"/> Chin Cup <i>Use HCPCS Code E2323 when Chin Cup not used in conjunction with a chin control.</i>	\$71	E2324
<input type="checkbox"/> Extended I-Handle	\$71	E2323
<input type="checkbox"/> Soft Ball Handle	\$50	E2323
<input type="checkbox"/> Medium Ball Handle	\$50	E2323
<input type="checkbox"/> Bodypoint 3 Inch Goal Post Handle	\$69	E2323
<input type="checkbox"/> Bodypoint 4 Inch Goal Post Handle	\$69	E2323


Joystick Mounting Brackets		
<input type="checkbox"/> Fixed, Inline <i>Not available with Q-Logic 2 EX Stand Alone Joystick.</i>	Standard	
<input type="checkbox"/> Left	<input type="checkbox"/> Right	
<input type="checkbox"/> Swing-Away	\$245	E1028
<input type="checkbox"/> Left	<input type="checkbox"/> Right	
<input type="checkbox"/> Inline	<input type="checkbox"/> Offset	
<input type="checkbox"/> Retract4 Swing-Away Joystick Arm Kit <i>Full 180° Movement</i>	\$250	E1028
<input type="checkbox"/> Left	<input type="checkbox"/> Right	
<input type="checkbox"/> Multi-Axis	\$100	E1028
<input type="checkbox"/> Left	<input type="checkbox"/> Right	
<input type="checkbox"/> No Joystick Mount		



Quantum Q6 Edge 2.0 with Contoured Seating

Description	Price	HCPCS Code
-------------	-------	------------

Programming

-  Q-Logic EX/NE Programmer \$545
- Q-Logic PC Programming Disk w/ Cable *Works with Q-Logic 2 EX only.* \$750

Contoured Seating

Depth-Adjustable High Back:

- Black Vinyl, Depth Adjustable High Back, Limited Recline No Charge
Seat comes automatically with a 60° belt and full length armrests. Any other selection would need to be added as an additional part.
- 18"W x 18-20"D 20"W x 18-20"D 22"W x 20-22"D

Limited Recline High Back:

- Black Vinyl, High Back with Limited Recline
- 16"W x 16"D 16"W x 20"D 18"W x 16"D 20"W x 16"D
- 16"W x 18"D 16"W x 22"D 18"W x 22"D 20"W x 22"D

Limited Recline High Back with Solid Seat Pan:

- Black Vinyl, High Back with Solid Seat Pan, Limited Recline No Charge
- 16"W x 16"D 16"W x 20"D 18"W x 16"D 18"W x 20"D 20"W x 16"D 20"W x 20"D
- 16"W x 18"D 16"W x 22"D 18"W x 18"D 18"W x 22"D 20"W x 18"D 20"W x 22"D

Medium Back:

- Black Vinyl, Medium Back No Charge
- Black Vinyl, Medium Back with Solid Seat Pan No Charge
- 16"W x 16"D 16"W x 18"D 18"W x 16"D 18"W x 18"D

Special Order Seat
 _____"W x _____"D
Seat size can not exceed 22"W x 22"D
Contact your inside sales specialist for pricing

Armrests

Armrest Type:

- | Description | Price | HCPCS Code |
|---|-------|------------|
| <input type="checkbox"/> Flip up, Height-Adjustable, Removable (each) | \$85 | E0973 |

Left Armpad:

- Full Length Desk Length¹

Right Armpad:

- Full Length Desk Length¹

1. Not available with depth-adjustable high back seat.

Power Elevation Seat/Towers

Description	Price	HCPCS Code
<input type="checkbox"/> Power Elevating Seat Actuator <i>SP code is required</i> <i>Not available with manual recliner. Max seat size is 22"W x 22"D. increases seat to floor height. Includes 1/4 speed inhibit.</i>	\$1,895	E2300
<input type="checkbox"/> Through Toggle, Left	Standard	
<input type="checkbox"/> Through Toggle, Right	Standard	
<input type="checkbox"/> Through Joystick	\$1,895	
Actuator:		
<input type="checkbox"/> 6" Stroke <input type="checkbox"/> 7" Stroke <i>Required on a 22"W contoured seat.</i>		
<input type="checkbox"/> Seat Towers		
<input type="checkbox"/> Standard Seat Mounts	Standard	
<input type="checkbox"/> Extended Seat Mounts	No Charge	

Contoured Seat to Floor Height

Contoured Seating:

- 20" Seat to Floor 21" Seat to Floor 22" Seat to Floor

Contoured Seating with Extended Towers:

- 23" Seat to Floor 24" Seat to Floor 25" Seat to Floor



Quantum Q6 Edge 2.0 with Contoured Seating

Description	Price	HCPCS Code
Lap Belts		

Description	Price	HCPCS Code
<input type="checkbox"/> 38" Lap Belt	Standard	
<input type="checkbox"/> 50" Lap Belt	Standard	
<input type="checkbox"/> 60" Lap Belt	Standard	
<input type="checkbox"/> 70" Lap Belt	Standard	
<input type="checkbox"/> 80" Lap Belt	Standard	
<input type="checkbox"/> 90" Lap Belt	Standard	
<input type="checkbox"/> 60" Padded Push Button, center Pull	\$115	
<input type="checkbox"/> Large 13-19", Hip Padded, Push Button	\$92	

Non-Powered Legrests		
-----------------------------	--	--

All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

<input type="checkbox"/> 60-90° Swing-Away Legrests adjustable in 5° Increments	\$162	
Rigging Selection:		
<input type="checkbox"/> Standard Matching Legrests (pair)	<input type="checkbox"/> Non-Matching Legrests (pair)	<input type="checkbox"/> Individual Legrests <i>Left or Right Only</i>
Footplate Selection:		
<input type="checkbox"/> Angle-Adjustable Footplates <i>Includes Heel Loops.</i>	\$75 ea.	K0040
<input type="checkbox"/> Standard, Extruded Footplates	No Charge	

Left Legrest:			Right Legrest:		
Lower Extensions:			Lower Extensions:		
<input type="checkbox"/> 11"-13" ¹	<input type="checkbox"/> 13"-16.5" ¹	<input type="checkbox"/> 15.5"-19.5" ¹	<input type="checkbox"/> 11"-13" ¹	<input type="checkbox"/> 13"-16.5" ¹	<input type="checkbox"/> 15.5"-19.5" ¹
<input type="checkbox"/> Angle Adjustable/Multi-Axis Footplate Size:			<input type="checkbox"/> Angle Adjustable/Multi-Axis Footplate Size:		
<input type="checkbox"/> 4"W x 6"D	<input type="checkbox"/> 5"W x 8"D	<input type="checkbox"/> 6"W x 11" ^{D3}	<input type="checkbox"/> 4"W x 6"D	<input type="checkbox"/> 5"W x 8"D	<input type="checkbox"/> 6"W x 11" ^{D3}
<input type="checkbox"/> 5"W x 6"D	<input type="checkbox"/> 6"W x 8" ^{D3}	<input type="checkbox"/> 8"W x 11" ^{D4}	<input type="checkbox"/> 5"W x 6"D	<input type="checkbox"/> 6"W x 8" ^{D3}	<input type="checkbox"/> 8"W x 11" ^{D4}
<input type="checkbox"/> Standard, Extruded Footplate Size:			<input type="checkbox"/> Standard, Extruded Footplate Size:		
<input type="checkbox"/> 7.25"W x 6"D	<input type="checkbox"/> 8.5"W x 6" ^{D3}		<input type="checkbox"/> 7.25"W x 6"D	<input type="checkbox"/> 8.5"W x 6" ^{D3}	

Accessories:			Price	HCPCS Code
<input type="checkbox"/> Heel Loops (pair) <i>Standard on Angle Adjustable and Multi-Axis footplates</i>			\$38	E0951
<input type="checkbox"/> Wheel Bumpers (pair)			\$25	
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair)			\$80	K0108
<input type="checkbox"/> Therafin Multi-Axis Amputee Support (each) <input type="checkbox"/> Left <input type="checkbox"/> Right <i>Swing-Away legrest upper assembly and receiver is required</i>			\$230	E1020
<input type="checkbox"/> Amputee Support Gel Cover (each) <input type="checkbox"/> Left <input type="checkbox"/> Right			\$150	K0108

- Not available with multi-axis footplates.
- Not available w/ Angle Adjustable Footplates
- May cause interference on seat widths less than 18"
- May cause interference on seat widths less than 20"

60-90° Swing-Away Legrest Compatibility Matrix	Contoured Seating		
	Low	Medium	High
Standard, 13" - 16.5"	!	✓	✓
Standard, 15.5" - 19.5"	!	!	!
Angle-Adjustable, 13" - 16.5"	!	✓	✓
Angle-Adjustable, 15.5" - 19.5"	!	!	!
Mul -Axis, 15.5" - 19.5"	!	!	!

- ✓ This rigging Configuration causes no interference issues at this seat-to-floor height.
- ! This rigging configuration may interfere with the ground (less than 1.75" clearance when fully extended) and will reduce obstacle climbing ability.



Quantum Q6 Edge 2.0 with Contoured Seating

Non-Powered Legrests (cont.)

All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

90° Pediatric Swing-Away Legrests No Charge

Rigging Selection:

Standard Matching Legrests (pair) Non-Matching Legrests (pair) Individual Legrests *Left or Right Only*

Footplate Selection:

Angle-Adjustable Footplates *Includes heel loops.* \$75 ea. K0040

Left Legrest:

Lower Extensions:

7"-10"L

Right Legrest:

Lower Extensions:

7"-10"L

Angle Adjustable Footplate Size:

4"W x 6"D 5"W x 6"D

4"W x 6"D

5"W x 6"D

Accessories:

Wheel Bumpers (pair) \$25

Pediatric Swing-Away Compatibility Matrix	Contoured Seating		
	Low	Medium	High
7"-10" Extension	✓	✓	✓

This rigging Configuration causes no interference issues at this seat-to-floor height.

Elevating Legrests \$325 E0990

Rigging Selection:

Standard Matching Legrests (pair) Non-Matching Legrests (pair) Individual Legrests *Left or Right Only*

Footplate Selection:

Angle-Adjustable Footplates *Includes Heel Loops.* \$75 ea. K0040

Standard, Extruded Footplates

Legrest Lower Extensions:

11"-14"L¹ 14"-17"L¹
 12"-14.5"L² 14.5"-18"L²

Right Legrest Lower Extensions:

11"-14"L¹ 14"-17"L¹
 12"-14.5"L² 14.5"-18"L²

Angle Adjustable Footplate Size:

Left Legrest:

4"W x 6"D 5"W x 8"D 6"W x 11"D³
 5"W x 6"D 6"W x 8"D³ 8"W x 11"D⁴

Right Legrest:

4"W x 6"D 5"W x 8"D 6"W x 11"D³
 5"W x 6"D 6"W x 8"D³ 8"W x 11"D⁴

Standard, Extruded Footplate Size:

Left Legrest:

7.25"W x 6"D 8.5"W x 6"D³

Left Legrest:

7.25"W x 6"D 8.5"W x 6"D³

Accessories:

Description

Heel Loops (pair) *Standard on Angle Adjustable and Multi-Axis footplates*

Price **HCPCS Code**

\$38 E0951

Wheel Bumpers (pair)

\$25

Gel Padded Adductor Bumpers (pair)

\$125 K0108

- 1. Not available with angle adjustable footplates.
- 2. Not available with standard, extruded footplates.

- 3. May cause interference on seat widths less than 18".
- 4. May cause interference on seat widths less than 20".

Elevating Legrest Compatibility Matrix	Contoured Seating		
	Low	Medium	High
Standard, 11" - 14"	✓	✓	✓
Standard, 14" - 17"	!	!	!
Angle-Adjustable, 12" - 14.5"	✓	✓	✓
Angle-Adjustable, 14.5" - 18"	!	!	!






This rigging Configuration causes no interference issues at this seat-to-floor height.

This rigging configuration may interfere with the ground (less than 1.75" clearance when fully extended) and will

This rigging configuration does not work with this seat-to-floor height.

Quantum Q6 Edge 2.0 with Contoured Seating

Description	Price	HCPCS Code
Additional Accessories		
<input type="checkbox"/> Oxygen Holder	\$150	E2208
<input type="checkbox"/> Weather Cover, Medium,Black	\$158	
<input type="checkbox"/> Cup Holder	\$17	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Rear Basket	\$103	
<input type="checkbox"/> Walker Holder	\$92	
<input type="checkbox"/> Cane/Crutch Holder	\$98	E2207
<input type="checkbox"/> Saddle Bag	\$28	
<input type="checkbox"/> Flag Holder	\$17	

Additional Foot Rigging Accessories		
<i>See Positioning Components order form for more options</i>		
<input type="checkbox"/> 	POSAlinc Knee Adductor for Swing-Away Legrests (pair)	\$512 K0108
<input type="checkbox"/> 	Stealth Foot Boxes (pair) <input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large	\$475 K0108
<input type="checkbox"/> 	Stealth Padded Shoe Holders (pair) <input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> X-Large	\$266 K0108
<input type="checkbox"/> 	Ankle Huggers (pair) <input type="checkbox"/> Medium <input type="checkbox"/> X-Large	\$156 K0108
<input type="checkbox"/> 	Bodypoint High Mount, Clamp-On Footplates (pair) <input type="checkbox"/> Angle Adjustable: <i>(Swing-away legrest upper assembly and receiver is required)</i> <input type="checkbox"/> Small, 4"x 7.75" <input type="checkbox"/> Large, 4.25" x 9.25"	\$356 K0040 \$356 K0040
<input type="checkbox"/>	<input type="checkbox"/> Multi-Axis <input type="checkbox"/> Large, 4.25" x 9.25"	\$460 K0040

Synergy® Cushions		
<input type="checkbox"/> Matching Seat Size	<input type="checkbox"/> Non-Matching Seat Size	
<i>Please refer to Synergy Cushion order forms for detailed pricing, sizes, codes and part #'s</i>		

Simplicity Cushion: General Use

<input type="checkbox"/> 16" x 16" - 20" x 24"	\$98	E2601
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$98	E2602

Solution® Cushion: Skin Protection & Positioning

<input type="checkbox"/> 10" x 10" - 20" x 24"	\$380	E2607
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$430	E2608

Solution® 1 Cushion: Skin Protection

<input type="checkbox"/> 14" x 14" - 20" x 24"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$430	E2604

Spectrum® Gel Cushion: Positioning

<input type="checkbox"/> 12" x 10" - 20" x 24"	\$480	E2605
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2606

Spectrum® Foam Cushion: Skin Protection & Positioning

<input type="checkbox"/> 10" x 10" - 20" x 24"	\$480	E2607
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2608

Spectrum® Air Flat Cushion: Adjustable Skin Protection

<input type="checkbox"/> 14" x 12" - 20" x 24"	\$480	E2622
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2623

Air Tube:

<input type="checkbox"/> Left Side	<input type="checkbox"/> Right Side
------------------------------------	-------------------------------------



Quantum Q6 Edge 2.0 with Contoured Seating

Description	Price	HCPCS Code
Synergy® Cushions (cont.)		
<i>Please refer to Synergy Cushion order forms for detailed pricing, sizes, codes and part #'s</i>		
Spectrum® Air Contour Cushion: Adjustable Skin Protection & Positioning		
<input type="checkbox"/> 14" x 12" - 20" x 24"	\$499	E2624
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$549	E2625
Air Tube:		
<input type="checkbox"/> Left Side		
<input type="checkbox"/> Right Side		
TRU-Comfort™2 4-Way Stretch: Skin Protection		
<input type="checkbox"/> 16" x 14" - 20" x 24"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 22" x 24"	\$480	E2604
TRU-Comfort™2 Air Mesh: Skin Protection		
<input type="checkbox"/> 16" x 14" - 20" x 24"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 22" x 24"	\$480	E2604
TRU-Comfort™ Accessory:		
<input type="checkbox"/> Sport Seat Cushion Cover	\$100	
<i>Available as an additional cover option for both TRU-Comfort seat cushion styles. Only available in sizes listed.</i>		
<input type="checkbox"/> 16"W x 16"D		
<input type="checkbox"/> 18"W x 18"D		
<input type="checkbox"/> 20"W x 20"D		