



Roscoe Mobility Quad Cane

The Roscoe Quad Cane is designed with a 4 point base, combined with an offset handle that provides additional stability and support for the patient.

Product Features

- 4 point base, combined with offset handle provides additional stability and support
- Easy-to-use, one-button height adjustment with locking anti-rattle ring allowing for quick, safe fitting
- Comfortable, contoured handle
- Attractive, easy-to-clean aluminum frame
- Available in attractive black finish
- Handle height adjusts from 29" to 38", in 1" increments
- Weight capacity: 300 lbs.
- Limited Lifetime Warranty
- HCPCS Code: E0105*

Item Number	Description	Color	Case Quantity
ROS-CN11050D-CH	Quad Cane, Small Base	Black	6
ROS-CN11070D-CH	Quad Cane, Large Base	Black	6

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com