



## *Roscoe Mobility* Quad Cane

The Roscoe Quad Cane is designed with a 4 point base, combined with an offset handle that provides additional stability and support for the patient.

### *Product* Features

- 4 point base, combined with offset handle provides additional stability and support
- Easy-to-use, one-button height adjustment with locking anti-rattle ring allowing for quick, safe fitting
- Comfortable, contoured handle
- Attractive, easy-to-clean aluminum frame
- Available in attractive black finish
- Handle height adjusts from 29" to 38", in 1" increments
- Weight capacity: 300 lbs.
- Limited Lifetime Warranty
- HCPCS Code: E0105\*

Item Number	Description	Color	Case Quantity
ROS-CN11050D-CH	Quad Cane, Small Base	Black	6
ROS-CN11070D-CH	Quad Cane, Large Base	Black	6

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



**To Order:**  
 (800) 3-ROSCOE  
 fax: (440) 572-4261  
 roscoemedical.com