



Roscoe Mobility Offset Aluminum Cane

The Roscoe Offset Aluminum Cane features a comfortable, contoured handle with a comfort grip, complete with convenient wrist strap. The attractive black finish resists marring and scratching.

Product Features

- Comfortable, contoured handle with comfort grip and strap
- Manufactured with sturdy 1" diameter anodized extruded aluminum tubing
- Handle height adjust in 1" increments
- Locking ring ensures added safety after adjustment
- Attractive finish resists scratched
- Handle height adjusts from 29" to 38"
- Weight capacity: 300 lbs.
- Limited Lifetime Warranty
- HCPCS Code: E0100*

Item Number	Description	Color	Case Quantity
ROS-CN10051B-10	Offset Cane with Strap	Black	10

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com